

Biochemical Diagnostics Detectabuse Liquid Urine Controls for 2015 SAMHSA CUTOFFS

- Human urine-based matrix simulates patient samples
- Up to 15 drug constituents within one control
- Stable for up to 3 years frozen, or refrigerated (unopened)
- Series XV-H in ready to use 20 mL vial size.
- Confirm 3 in ready to use 50 mL vial size.
- Each Series XV-H control contains: THC, Benzoylecgonine, PCP, Morphine, 6-MAM, d-Methamphetamine, MDMA, Secobarbital, Oxazepam, Methadone, Methaqualone, Propoxyphene, Oxycodone and Hydrocodone
- Each Confirm 3 control contains: THC, Benzoylecgonine PCP, Codeine, Morphine (OPI), 6-MAM, d-Amphetamine d-Methamphetamine, MDMA, MDA, MDEA, Hydrocodone, Hydromorphone, Oxycodone, Oxymorphone



The DETECTABUSE® Advantage:

Detectabuse® Series XV-H & Confirm 3 controls are designed to monitor the performance of qualitative or quantitative procedures for the detection of drugs in urine. We also offer a complete line of controls for other immunoassay, DOA devices, GC/MS and custom formulations.

NEW! DETECTABUSE® LIQUID CONTROL URINE IMMUNOASSAY SERIES XV-H AND GC/MS CONFIRM 3		
Catalog #	Series XV-H (High Opiate, 2000 ng/mL Cutoff)	Vial Size
19601057	Series XV-H Cutoff -25%	20 mL
19601058	Series XV-H Cutoff +25%	20 mL
19601059	Series XV-H 3X Cutoff *special order	20 mL
GC/MS Confirm 3 (High Opiate, 2000 ng/mL cutoff)		
19582040	Confirm 3 Cutoff -60%	50 mL
19582025	Confirm 3 Cutoff -25% *special order	50 mL
19582000	Confirm 3 Cutoff	50 mL
19582125	Confirm 3 Cutoff +25%	50 mL



Contact: Colleen Gang, Product Manager, Liquid Controls • 180 Heartland Blvd, Edgewood, NY 11717
• E-mail: cgang@biochemicaldiagnostics.com • Phone: 1-800-223-4835 • Fax: 1-631-595-9204