

KOVA INTERNATIONAL



KOVA-Trol™ and KOVA Liqua-Trol™ Shelf Life Extension for All New Lots

Key Message

Kova International has announced extended shelf life stability for its wide range of products. These products include KOVA-Trol I™ with and without UB and KOVA-Trol II™. KOVA® Liqua-Trol™ Normal and Abnormal Urinalysis controls (with and without microscopic). The shelf life stability has been increased from 18 months to 27 months from the date of manufacture. Clinical laboratory customers can now rest assured that the KOVA® Urinalysis Controls with extended stability will meet and exceed their expectations for quality assurance of their Urinalysis testing requirements.

Call to Action

Kova International understands the critical importance of accuracy and dependability of Urinalysis supplies. Extensive stability and validation testing has been completed to insure that the longer timeframes will have no impact on the quality or dependability of products. As a result of the increased shelf life, testing sites may experience a reduction in lot turnover and the associated costs.

KOVA-Trol™ is still being offered in three levels to support monitoring across the entire decision range of reagent strip chemistries. The three levels are High Abnormal, which has optional Urobilinogen values, Low Abnormal, which contains creatinine, and microalbumin, and Normal, which contains hCG along with creatinine, and Microalbumin.

Features and Benefits:

- With proper storage and handling, our traditional KOVA-Trol I™, and KOVA-Trol II™ and KOVA-Trol III™ remain in a stable state for up to 27 months from the date of manufacture.
- When KOVA-Trol I™, KOVA-Trol II™, and KOVA-Trol III™ are reconstituted, all constituents are stable from day 1 of reconstitution until day 7 at 2-8° C.
- KOVA® Liqua-Trol™ Normal and Abnormal, are ready to use liquid for urine control testing. It is available both with and without microscopics.

KOVA® A Complete Urinalysis System



The KOVA® System includes specially engineered microscopic slides, tubes & caps, custom pipettes or “petters”, collection cups and multi-level urinalysis controls for reliable, simple and effective workflow standardization. From specimen collection to slide preparation, KOVA® Urinalysis products help ensure samples are properly collected and handled for chemical tests using urine dipstick chemistry, confirmatory and pregnancy tests, and microscopic urinalysis. They give you the control you need for reliable, reproducible results with minimal error and maximum efficiency.

- Disposable components and quality control materials.
- Convenient disposable microscope slides with optional counting grid.
- Multilevel urine chemistry and microscopy controls in convenient lyophilized and liquid formats.
- Multilevel serum protein controls for refractometry.

Resources

Samples: please email cseurope@kovaintl.com with subject line: KOVA SAMPLE PACK

Website: www.kovaintl.com

Sales/Technical Support: +31 (0) 20 658 6162